



Electrodes

Valid from 1st July, 2023
All Prices are in Australian Dollars and do NOT include GST or freight.

Prices are subject to change without notification.

Ionode Pty Ltd

12 Walker Street
Tennyson QLD 4105 Australia
Tel: +61 7 3848-1660
Email: sales@ionode.com
Web: www.ionode.com

EFT Payment Details

Bank: Westpac Banking Corp.
Account Name: Ionode Pty Ltd
BSB: 034-037
Account Number: 252725
SWIFT Code: WPACAU2S

Product	Order Code	Description	Price	Extra Cable \$/m
pH Electrodes	Ionode's unique intermediate junction provides best results in almost any sample including food & dairy, cosmetics, wine & mining. Submersible & suitable for in-line use with MK01 kit or similar. Standard sensor includes 1m cable, BNC connector, 2 junction sleeves & 45ml reference electrolyte.			
IJ40	1111A-01.0-B	pH sensor with bullet membrane for general use	\$ 299.00	\$ 10.00
IJ40-ALK	1111C-01.0-B	IJ40 sensor suitable for low sodium error applications (for pH>12 in presence of Na ions)	\$ 309.00	\$ 10.00
IJ40-HT	1211A-01.0-B	IJ40 sensor suitable for high temp (up to 100°C). Solvent resistant	\$ 389.00	\$ 10.00
IJ40 ATC	8111A-01.0-BXX-X	IJ40 with built in ATC - contact us for quote/specifications	\$ 399.00	\$ 17.50
IJ40 HT ATC	8211A-01.0-BXX-X	IJ40 HT with built in ATC - contact us for quote/specifications	\$ 489.00	\$ 17.50
IJ44 HT ATC	8215A-01.0-BXX-X	IJ44 HT with built in ATC - contact us for quote/specifications	\$ 489.00	\$ 17.50
IJ44	1115A-01.0-B	pH sensor with spear membrane ideal for food/dairy applications, creams, emulsions etc.	\$ 299.00	\$ 10.00
IJ44-HT	1215A-01.0-B	IJ44 sensor suitable for high temp (up to 100°C). Solvent resistant	\$ 389.00	\$ 10.00
IJ46	1114A-01.0-B	pH sensor with flat membrane suitable for flat surface, fabric measurement etc.	\$ 389.00	\$ 10.00
IJ46AT	8114A-01.0-BXX-X	IJ46 with built in ATC - contact us for quote/specifications	\$ 489.00	\$ 17.50
pH Electrodes	Double junction, sealed gel electrodes for routine applications with low contaminant levels and reasonable ionic strength. Submersible and suitable for in-line use with MK01 kit or similar. Standard supply is 1m cable, BNC connector.			
IH40	1131A-01.0-B	Sealed gel pH sensor with bullet membrane for general use	\$ 173.00	\$ 10.00
IH40 ATC	8131A-01.0-BXX-X	IH40 with built in ATC - contact us for quote/specifications	\$ 273.00	\$ 17.50
IH40-ALK	1131C-01.0-B	IH40 suitable for low sodium error applications (for pH>12 in presence of Na ions)	\$ 183.00	\$ 10.00
pH Electrodes	P Series - Single junction, refillable electrodes for routine applications, education. Standard supply is 1m cable, BNC connector.			
PBFA	1142A-01.0-B	Single junction (Ag/AgCl) refillable pH sensor for general purpose	\$ 150.00	\$ 10.00
PBFA ATC	8142A-01.0-BXX-X	PBFA with built in ATC - contact us for quote/specifications	\$ 250.00	\$ 17.50
PBFC	1142C-01.0-B	PBFA sensor for low sodium error applications (for pH>12 in presence of Na ions)	\$ 160.00	\$ 10.00
Reference Electrodes	Various junction types. Standard supply is 1m cable, 4mm connector.			
IJ14	4119-01.0-4	Intermediate junction reference electrode	\$ 210.00	\$ 5.00
IH10	4139-01.0-4	Sealed gel double junction reference electrode	\$ 130.00	\$ 5.00
PJFO	4149-01.0-4	Refillable single junction reference electrode	\$ 135.00	\$ 5.00
Silver/Silver Chloride	GBDJ-AgAgCl	Glass-bodied double junction Ag/AgCl reference electrode	\$ 263.00	\$ 5.00
Silver/Silver Sulfate	GBDJ-AgAgSu	Glass-bodied double junction Ag/Ag2SO4 reference electrode	\$ 315.00	\$ 5.00
Silver/Silver Oxide - NEW	PBDJ-AgAg2O	PTCFE double junction Ag/Ag2O reference electrode	\$ 368.00	\$ 5.00
Calomel	Cal-Ref-01-K	Glass-bodied refillable calomel reference electrode with porous frit	\$ 295.00	\$ -



Electrodes

Valid from 1st July, 2023
All Prices are in Australian Dollars and do NOT include GST or freight.

Prices are subject to change without notification.

Ionode Pty Ltd

12 Walker Street
Tennyson QLD 4105 Australia
Tel: +61 7 3848-1660
Email: sales@ionode.com
Web: www.ionode.com

EFT Payment Details

Bank: Westpac Banking Corp.
Account Name: Ionode Pty Ltd
BSB: 034-037
Account Number: 252725
SWIFT Code: WPACAU2S

Product	Order Code	Description	Price	Extra Cable \$/m
ORP/Metal Electrodes	Featuring Ionode's unique intermediate junction assuring reliable results in just about any sample type including mining and electroplating. Submersible & suitable for in-line use with MK01 kit or similar. Standard supply is 1m cable, BNC connector, 2 reference sleeves & 45ml bottle of reference electrolyte.			
IJ64-Pt	211Pt-01.0-B	Intermediate Junction, Pt wire ORP. All applications ranging from low ionic to harsh or highly contaminated	\$ 320.00	\$ 6.00
IJ64-Pt-Disk	211PtD-01.0-B	Pt disc - for faster response in various samples	\$ 335.00	\$ 6.00
IJ64-Au	211Au-01.0-B	Intermediate Junction, Au wire ORP. All applications ranging from low ionic to harsh or highly contaminated	\$ 351.00	\$ 6.00
IJ64-HT	221Pt-01.0-B	Pt wire ORP - suitable for high temp (to 100°C). Solvent resistant	\$ 410.00	\$ 6.00
IJ-Ag	211Ag-01.0-B	Silver billet electrode for silver nitrate titrations for salt etc	\$ 331.00	\$ 6.00
ORP/Metal Electrodes	Double junction, sealed gel electrodes for routine applications with low contaminant levels and reasonable ionic strength. Submersible and suitable for in-line use with MK01 kit or similar. Standard supply is 1m cable, BNC connector.			
IH30-Pt	213Pt-01.0-B	Pt wire ORP - double junction sealed gel	\$ 195.00	\$ 6.00
IH30-Pt-Disk	213PtD-01.0-B	Pt disc - for faster response	\$ 210.00	\$ 6.00
IH30-Au	213AuD-01.0-B	Au ORP sensor	\$ 230.00	\$ 6.00
ORP/Metal Electrodes	P Series - Single junction, refillable electrodes for routine applications, education. Standard supply is 1m cable with BNC connector.			
PRFO	4149-01.0-B	Pt wire ORP	\$ 135.00	\$ 6.00
ISE Electrodes	Refillable ion selective electrodes, combination type (double junction) with epoxy body, ceramic frit, 1m cable with BNC connection.			
Ammonia (gas sensing)	332NH3-01.0-B	Ammonia ion selective electrode with 2 spare membrane caps	\$ 780.00	\$ 10.00
Bromide	364Br-01.0-B	Bromide ion selective electrode (Blue Series). Improved detection limit & quicker response	\$ 780.00	\$ 10.00
Chloride	364Cl-01.0-B	Chloride ion selective electrode (Blue Series). Improved detection limit & quicker response	\$ 780.00	\$ 10.00
Cyanide	364CN-01.0-B	Cyanide ion selective electrode (Blue Series). Improved detection limit & quicker response	\$ 780.00	\$ 10.00
Fluoride	364F-01.0-B	Fluoride ion selective electrode (Blue Series). Improved detection limit & quicker response	\$ 830.00	\$ 10.00
Iodide	364I-01.0-B	Iodide ion selective electrode (Blue Series). Improved detection limit & quicker response	\$ 780.00	\$ 10.00
Nitrate	364NO3-01.0-B	Nitrate ion selective electrode (Blue Series). Improved detection limit & quicker response	\$ 780.00	\$ 10.00
Silver/silver sulphide	364Ag2S-01.0-B	Silver/silver sulphide ion selective electrode (Blue Series). Improved detection limit & quicker response	\$ 780.00	\$ 10.00
Calcium (PVC)	364Ca-01.0-B	Calcium ISE (Blue Series). Improved detection limit & quicker response with regenerable PVC membrane	\$ 780.00	\$ 10.00
Calcium Hardness (PVC)	364CH-01.0-B	Hardness ISE (Blue Series). Improved detection limit & quicker response with regenerable PVC membrane	\$ 780.00	\$ 10.00
Potassium (PVC)	364K-01.0-B	Potassium ISE (Blue Series). Improved detection limit & quicker response with regenerable PVC membrane	\$ 780.00	\$ 10.00
Sodium (PVC)	364Na-01.0-B	Sodium ISE (Blue Series). Improved detection limit & quicker response with regenerable PVC membrane	\$ 780.00	\$ 10.00
Surfactant (PVC)	364SF-01.0-B	Surfactant ISE (Blue Series). Improved detection limit & quicker response with regenerable PVC membrane	\$ 780.00	\$ 10.00
PVC Membrane Regeneration Service	PVC-ISE-ReGen	Only applicable for PVC membrane	\$ 130.00	\$ -



Electrodes

Valid from 1st July, 2023
All Prices are in Australian Dollars and do NOT include GST or freight.

Prices are subject to change without notification.

Ionode Pty Ltd

12 Walker Street
Tennyson QLD 4105 Australia
Tel: +61 7 3848-1660
Email: sales@ionode.com
Web: www.ionode.com

EFT Payment Details

Bank: Westpac Banking Corp.
Account Name: Ionode Pty Ltd
BSB: 034-037
Account Number: 252725
SWIFT Code: WPACAU2S

Product	Order Code	Description	Price	Extra Cable \$/m
ISE Electrodes	Refillable ion selective electrodes, combination type (double junction) with epoxy body, ceramic frit, 1m cable with BNC connection.			
IJ Bromide ISE	311Br-01.0-B	Intermediate junction bromide ion selective electrode - suitable for titration applications	\$ 899.00	\$ 10.00
IJ Chloride ISE	311Cl-01.0-B	Intermediate junction chloride ion selective electrode - suitable for titration applications	\$ 899.00	\$ 10.00
IJ Cyanide ISE	311CN-01.0-B	Intermediate junction cyanide ion selective electrode - suitable for titration applications	\$ 899.00	\$ 10.00
IJ Fluoride ISE	311F-01.0-B	Intermediate junction fluoride ion selective electrode - suitable for titration applications	\$ 949.00	\$ 10.00
IJ Iodide ISE	311I-01.0-B	Intermediate junction iodide ion selective electrode - suitable for titration applications	\$ 899.00	\$ 10.00
ISE Polishing Kit	ISE-PolKit	Kit with polishing cloth, glass backing plate and polishing solution	\$ 120.00	\$ -
Conductivity Probes	Conductivity Sensors, 2-pole. Standard supply is 1m cable, BNC connector. ATC versions BNC/RCA to suit Elmetron conductivity meters.			
ECF-1	ELECF-1	PVC body, K=0.45, for most samples, wide Range 0uS/cm - 500 mS/cm	\$ 369.00	\$ -
ECF-1T	ELECF-1T	PVC body, K=0.45, for most samples, wide Range 0uS/cm - 500 mS/cm with ATC	\$ 455.00	\$ -
Voltammetry Electrodes	Working electrodes (WE), Pt/Aux electrodes (CE), and Reference electrodes (RE) for electrochemistry studies. With pin connector.			
Pt WE 1mm	Echem-WE-Pt-1	PEEK body 3mm diameter with 1mm diameter Pt disc	\$ 252.00	\$ -
Au WE 1mm	Echem-WE-Au-1	PEEK body 3mm diameter with 1mm diameter Au disc	\$ 272.00	\$ -
GC WE 1mm	Echem-WE-GC-1	PEEK body 3mm diameter with 1mm diameter GC disc (glassy carbon)	\$ 282.00	\$ -
Pt WE 3mm	Echem-WE-Pt	PEEK body 6mm diameter with 3mm diameter Pt disc	\$ 252.00	\$ -
Au WE 3mm	Echem-WE-Au	PEEK body 6mm diameter with 3mm diameter Au disc	\$ 272.00	\$ -
GC WE 3mm	Echem-WE-GC	PEEK body 6mm diameter with 3mm diameter GC disc (glassy carbon)	\$ 282.00	\$ -
CE/Aux electrode	Echem-Aux-Pt	Glass body 6mm diameter counter electrode with 6x6mm Pt foil	\$ 252.00	\$ -
Ref-refillable	Echem-Ref-R	Refillable Ag/AgCl reference electrode. Supplied with small bottle of electrolyte	\$ 105.00	\$ -
Ref-refillable	Echem-Ref-R-SR	Refillable reference electrode. For non-aqueous measurements	\$ 105.00	\$ -
Microelectrode Carbon Fibre 7 μ	Echem-MWE-CF-7M	Glass bodied microelectrode for electrochemical research - contact us for custom options	\$ 662.00	\$ -
Microelectrode Platinum 10 μ	Echem-MWE-Pt-10	Glass bodied microelectrode for electrochemical research - contact us for custom options	\$ 515.00	\$ -
Microelectrode Platinum 25 μ	Echem-MWE-Pt-25	Glass bodied microelectrode for electrochemical research - contact us for custom options	\$ 515.00	\$ -
Additional Junction	Echem-Ref-A-J	Additional junction for EChem Reference	\$ 126.00	\$ -
EChem Polishing Kit	Echem-PolKit	Kit with polishing cloths, glass backing plates and polishing solutions	\$ 290.00	\$ -



Electrodes

Valid from 1st July, 2023
All Prices are in Australian Dollars and do NOT include GST or freight.

Prices are subject to change without notification.

Ionode Pty Ltd

12 Walker Street
Tennyson QLD 4105 Australia
Tel: +61 7 3848-1660
Email: sales@ionode.com
Web: www.ionode.com

EFT Payment Details

Bank: Westpac Banking Corp.
Account Name: Ionode Pty Ltd
BSB: 034-037
Account Number: 252725
SWIFT Code: WPACAU2S

Product	Order Code	Description	Price	Extra Cable \$/m
Process Light/Industrial Sensors	Industrial electrodes with detachable cable. Standard supply is S8/PG13.5 connector. (Cable Assemblies located on page 5)			
IJ40-S8	1111A-00.0-S8	Intermediate Junction pH sensor with bullet membrane for general use with S8/PG13.5 connector	\$ 329.00	\$ -
IJ64-Pt-S8	211Pt-00.0-S8	Intermediate Junction, Pt wire ORP with S8/PG13.5 connector. All applications ranging from low ionic to harsh or highly contaminated	\$ 350.00	\$ -
IH40-S8	1131A-00.0-S8	Double Junction, sealed gel pH sensor with bullet membrane for general use with S8/PG13.5 connector	\$ 203.00	\$ -
IH30-Pt-S8	213Pt-00.0-S8	Pt wire ORP - double junction sealed gel with S8/PG13.5 connector	\$ 225.00	\$ -
IH30-Pt-Disk-S8	213PtD-00.0-S8	Pt disc - for faster response - double junction sealed gel with S8/PG13.5 connector	\$ 240.00	\$ -
IHK1	5309-S8	K=1 Conductivity, 2-pole, abraded titanium pins. S8/PG13.5 connector	\$ 205.00	\$ -
Process Light/Industrial Sensors	In-line single junction industrial grade sensors. Standard supply is 1.8m cable, BNC connector.			
SI20	1322A-01.8-B	1/4" BSPM, single junction pH sensor for clean samples; pools, cooling towers, drinking water etc.	\$ 150.00	\$ 10.00
SI30	232Pt-01.8-B	1/4" BSPM, single junction Pt ORP sensor for clean samples; pools, cooling towers, drinking water etc.	\$ 163.00	\$ 6.00
CS20	1421A-01.8-B	3/4" BSPF, single junction pH sensor for clean samples; pools, cooling towers, drinking water etc.	\$ 168.00	\$ 10.00
CS30	242Pt-01.8-B	3/4" BSPF, single junction Pt ORP sensor for clean samples; pools, cooling towers, drinking water etc.	\$ 184.00	\$ 6.00
CS30-Au	242Au-01.8-B	3/4" BSPF, single junction Au ORP sensor for clean samples; pools, cooling towers, drinking water etc.	\$ 204.00	\$ 6.00
Multi Parameter Options	Contact Us	pH/ORP, temp, solution earth, conductivity options available.	Contact us	
Process Light/Industrial Sensor Accessories	Process/Light Industrial Sensors Accessories			
FT01	IN-FT01	3/4" x 3/4" BSPM Tee-piece adapter for CS series	\$ 7.50	\$ -
PS01	IN-PS01	Pool Saddle, 1/4" BSP. Only suitable for 40mm pipe	\$ 7.50	\$ -
Electrode Accessories	Electrode adapters, mounting accessories, stirrer stands & electrode holders.			
EH-G0001	GEH-G0001	Electrode Holder	\$ 165.00	\$ -
SS-G0001	GSS-G0001A	Magnetic Stirrer Stand	\$ 276.00	\$ -
BNC/DIN Adaptor	CABD	BNC Cable Adaptor to suit WTW and many German meters	\$ 74.00	\$ -
BNC/LEMO Adaptor	CABL	BNC Adaptor to suit Metrohm (5mm Lemo)	\$ 74.00	\$ -
MK01	MK01	Mounting Kit for IH/IJ Ionode Electrodes	\$ 58.00	\$ -



Electrodes

Valid from 1st July, 2023
All Prices are in Australian Dollars and do NOT include GST or freight.

Prices are subject to change without notification.

Ionode Pty Ltd

12 Walker Street
Tennyson QLD 4105 Australia
Tel: +61 7 3848-1660
Email: sales@ionode.com
Web: www.ionode.com

EFT Payment Details

Bank: Westpac Banking Corp.
Account Name: Ionode Pty Ltd
BSB: 034-037
Account Number: 252725
SWIFT Code: WPACAU2S

Product	Order Code	Description	Price	Extra Cable \$/m
Cable Assemblies		Cable Assemblies for S7/S8 Connector Electrodes		
S7/S8 3m Cable Assembly	22.030	S7/S8 Female Cable Assembly, 3m	\$ 93.50	\$ -
S7/S8 3m Cable Assembly - BNC	22.031 - BNC	S7/S8 Female Cable Assembly, 3m, BNC connector	\$ 115.00	\$ -
S7/S8 10m Cable Assembly	22.034	S7/S8 Female Cable Assembly, 10m	\$ 181.50	\$ -
S7/S8 10m Cable Assembly - BNC	22.034 - BNC	S7/S8 Female Cable Assembly, 10m, BNC connector	\$ 203.00	\$ -
S7/S8 Custom Cable Assembly - BNC	22.03x - BNC	S7/S8 custom cable assembly/length - 1m or 2m BNC	\$ 195.00	\$ -
S7/S8 Custom Cable Assembly - BNC	22.03x - BNC	S7/S8 custom cable assembly/length - 4m to 9m BNC	\$ 195.00	\$ -
Electrode Upgrades		Additional upgrades for electrodes		
Temperature Upgrade (pH)	Temp_pH	Addition of PT100, PT1000, 10K or 30K temperature sensor to pH electrode	\$ 100.00	\$ -
Temperature Upgrade (ORP)	Temp_ORP	Addition of PT100, PT1000, 10K or 30K temperature sensor to ORP electrode	\$ 80.00	\$ -
Platinum Pin to Disk Upgrade	Disk	Upgrade platinum pin to disk in ORP electrode	\$ 15.00	\$ -
Precious Metal - Gold Upgrade	Au	Upgrade platinum to gold in ORP electrode	\$ 20.00	\$ -
High Temp Material Upgrade	HT	High temp upgrade (up to 100°C). Solvent resistant	\$ 90.00	\$ -
Connector Upgrades		Rework of standard BNC to another connector.		
Upgrade BNC to LEMO	LEMO	Rework BNC and replace with LEMO	\$ 35.00	\$ -
Upgrade BNC to DIN	DIN	Rework BNC and replace with DIN	\$ 35.00	\$ -
Upgrade BNC to S8	S8	Rework BNC and replace with S8	\$ 30.00	\$ -
Upgrade BNC to Crimp Pins	PINS	Rework BNC and replace with Crimp Pins	\$ 10.00	\$ -



Valid from 1st July, 2023
All Prices are in Australian Dollars and do NOT include GST or freight.

Prices are subject to change without notification.

12 Walker Street
Tennyson QLD 4105 Australia
Tel: +61 7 3848-1660
Email: sales@ionode.com
Web: www.ionode.com

Bank: Westpac Banking Corp.
Account Name: Ionode Pty Ltd
BSB: 034-037
Account Number: 252725
SWIFT Code: WPACAU2S

Product	Order Code	Description	Price
pH Buffers		pH buffers, with dye for easy identification	
pH 4.01/250ml	pH4-250	pH 4.01 buffer 250mL, red	\$ 19.30
pH 4.01/ 1L	pH4-1L	pH 4.01 buffer 1L, red	\$ 48.20
pH 6.86 / 250mL	pH6.86-250	pH 6.86 buffer 250mL, purple	\$ 24.00
pH 6.86 / 1L	pH6.86-1L	pH 6.86 buffer 1L, purple	\$ 48.20
pH 7.00 / 250mL	pH7-250	pH 7.00 buffer 250mL, green	\$ 18.10
pH 7.00 / 1L	pH7-1L	pH 7.00 buffer 1L, green	\$ 48.20
pH 10.01 / 250mL	pH10-250	pH 10.01 buffer 250 mL, blue	\$ 20.30
pH 10.01 / 1L	pH10-1L	pH 10.01 buffer 1L, blue	\$ 48.20
Conductivity Standards		KCl standard solutions for conductivity calibration. Other values available on request.	
1,413 µS/cm, 250ml	EC1413-250	Conductivity standard, 250mL	\$ 50.30
1,413 uS/cm, 1L	EC1413-1L	Conductivity standard, 1L	\$ 78.10
2,760 µS/cm, 1L	EC2760-1L	Conductivity standard, 1L	\$ 78.10
12,880 µS/cm, 500ml	EC12880-500	Conductivity standard, 500mL	\$ 64.20
ORP Standards		ORP standard solutions for checking ORP electrodes	
Zobells ORP Std "A" & "B" 250ml	ORP220-250	ORP standard, Part A & Part B, 229mV @ 25°C	\$ 160.50
EasyCal ORP 700, 250ml	ORP700-250	ORP 700mV EasyCal Standard, 250mL	\$ 35.30
EasyCal ORP 700, 1L	ORP700-1L	ORP 700mV EasyCal Standard, 1L	\$ 114.50
Reference Gels & Electrolytes		Reference electrolyte solutions and gels	
RE45	RE45	KCl reference electrolyte gel, 45mL in squeeze bottle - to suit IJ series electrodes	\$ 49.00
KN45	KN45	KNO3 reference electrolyte, 45mL in squeeze bottle	\$ 32.00
KN45G	KN45G	KNO3 reference gel electrolyte, 45mL in squeeze bottle - to suit IJ series ISE's	\$ 53.00
RE30	RE30	3M KCl reference electrolyte gel, 30mL in squeeze bottle	\$ 32.00
3M KCl	KC45	3M KCl solution, 45mL in squeeze bottle	\$ 32.00
3M KCl.Ag	KC45S	3M KCl solution sat'd with AgCl, 45mL in squeeze bottle	\$ 32.00
Maintenance Solutions		Electrode storage solutions	
KS250	KS250	pH electrode storage solution, 250mL	\$ 53.00
KS1000	KS1000	pH electrode storage solution, 1L	\$ 126.00



Pocket/Handheld Meters

Valid from 1st July, 2023
All Prices are in Australian Dollars and do NOT include GST or freight.

Prices are subject to change without notification.

Ionode Pty Ltd

12 Walker Street
Tennyson QLD 4105 Australia
Tel: +61 7 3848-1660
Email: sales@ionode.com
Web: www.ionode.com

EFT Payment Details

Bank: Westpac Banking Corp.
Account Name: Ionode Pty Ltd
BSB: 034-037
Account Number: 252725
SWIFT Code: WPACAU2S

Product	Order Code	Description	Price
Pocket pH Meter for use in meats, dairy, semi-solids, creams, emulsions			
CPx-105	ELCPx-105A	Pocket pH/mV meter with IJ44 ATC sensor, RE30 gel in carry case	\$ 492.20
GPx-105	ELGPx-105	Replacement IJ44 pH/ATC head to suit CPx-105	\$ 363.80
Handheld Meters - Basic Models (in carry case with temp sensor, battery powered)			
CP-411	ELCP-411d	pH/mV meter, 3 point autocal with temp sensor	\$ 599.20
CP-411-IH Package	ELCP-411-IH	pH/mV meter, 3 point autocal with temp sensor and IH40A sealed gel pH sensor	\$ 733.60
CP-411-IJ Package	ELCP-411-IJ	pH/mV meter, 3 point autocal with temp sensor and IJ40A intermediate junction pH sensor	\$ 853.30
CPC-411 Package	ELCPC-411-Pkg	pH/EC meter with IH40A and ECF-1T	\$ 1,165.80
CC-411 Package	ELCC-411-Pkg	Conductivity meter with ECF-1T	\$ 1,063.10
CO-411 Package	ELCO-411-Pkg	DO meter with sensor - galvanic type, with 2 spare membranes and electrolyte	\$ 889.50
Handheld Meters - GLP Compliant - in carry case with power adapter, temp. sensor, & with data capture software and USB cable			
CP-401	ELCP-401d	CP-401 pH/mV handheld meter	\$ 840.00
CP-401-IH Package	ELCP-401-IH	CP-401 pH/mV handheld meter w/ IH40 sensor	\$ 1,003.00
CC-401	ELCC-401d	CC-401 EC handheld meter	\$ 1,171.70
CC-401 Package	ELCC-401-Pkg	CC-401 EC handheld meter with ECF-1T	\$ 1,545.40
CPC-401	ELCPC-401d	CPC-401 pH/EC handheld meter	\$ 1,385.70
CPC-401 Package	ELCPC-401-pkg	CPC-401 pH/EC handheld meter with IJ40A & ECF-1T	\$ 2,032.70
CX-401	ELCX-401d	CX-401 multi-function handheld meter	\$ 1,219.80
CX-401 Package	ELCX-401-Pkg	CX-401 multi-function handheld meter w/ ECF-1 and IJ pH sensor	\$ 1,823.90
Accessories			
COG-1 kit	ELCOG-1Kit	Replacement membrane kit for COG-1 (2 membranes, 30mL bottle of electrolyte)	\$ 130.00
COG-1	ELCOG-1	COG-1 DO sensor, galvanic to suit Elmetron DO meters	\$ 631.30
12V Power Adapter	30.009	SMP500-12R Power Adapter for 411 Series Meters	\$ 42.80
6V Power Adapter	30.005	SMP800-6R Power Adapter for 401 Series Meters	\$ 42.80
9V Power Adapter	30.003	SMP1200-9R Power Adapter for CSP-501 and 600 Series Meters	\$ 42.80
USB-charger	30.008	USB wall adapter/charger - excl USB cable	\$ 32.10
Temperature sensor	ELTempSensor	Temperature Sensor for Elmetron Meters (Pt1000)	\$ 176.60



Lab Meters

Valid from 1st July, 2023
All Prices are in Australian Dollars and do NOT include GST or freight.

Prices are subject to change without notification.

Ionode Pty Ltd

12 Walker Street
Tennyson QLD 4105 Australia
Tel: +61 7 3848-1660
Email: sales@ionode.com
Web: www.ionode.com

EFT Payment Details

Bank: Westpac Banking Corp.
Account Name: Ionode Pty Ltd
BSB: 034-037
Account Number: 252725
SWIFT Code: WPACAU2S

Product	Order Code	Description	Price
Bench Top - Basic models with temp sensor and power adapter			
CP-511	ELCP-511d	CP-511 pH/mV bench top meter	\$ 1,066.10
CP-511 Package	ELCP-511-Pkg	CP-511 bench top meter with IJ40A sensor	\$ 1,392.00
CC-511	ELCC-511d	CC-511 EC bench top meter	\$ 1,230.50
CC-511 Package	ELCC-511-Pkg	CC-511 bench top meter with ECF-1	\$ 1,601.20
CPC-511	ELCPC-511d	CPC-511 pH/EC bench top meter	\$ 1,524.80
CPC-511 Package	ELCPC-511-Pkg	CPC-511 bench top meter with IJ40A sensor and ECF-1	\$ 2,083.20
Bench Top Meters - Advanced models (GLP compliant) - with power adapter, temperature sensor & data collection software			
CPI-505	ELCPI-505d	CPI-505 pH/Ion bench top meter - without sensor	\$ 1,599.70
CX-601	ELCX-601d	CX-601 multi function bench top meter - without sensor	\$ 2,113.30